

## Product Datasheet

### SLC25A4 antibody (orb681303)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | SLC25A4 antibody  |
| <b>Species/Host</b>        | Rabbit  |
| <b>Reactivity</b>          | Human, Mouse, Rat   |
| <b>Conjugation</b>         | Unconjugated  |
| <b>Tested Applications</b> | ELISA, WB   |
| <b>Immunogen</b>           | Human SLC25A4   |
| <b>Target</b>              | SLC25A4   |
| <b>Preservatives</b>       | PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles. |
| <b>Form/Appearance</b>     | Liquid  |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                         |
| <b>Note</b>                | For research use only   |
| <b>Isotype</b>             | IgG   |
| <b>Clonality</b>           | Polyclonal  |
| <b>Uniprot ID</b>          | <a href="#">P12235</a>  |
| <b>Expiration Date</b>     | 12 months from date of receipt.   |