

Product Datasheet

HSPG2 antibody (orb682517)

| | |
|----------------------------|---|
| Description | HSPG2 antibody |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, WB |
| Immunogen | Human HSPG2-Specific |
| Target | HSPG2 |
| Preservatives | PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles. |
| Form/Appearance | Liquid |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | P98160 |
| Expiration Date | 12 months from date of receipt. |