


Product Datasheet

FIGF antibody (orb682864)

| | |
|----------------------------|--|
| Description | FIGF antibody |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, IHC, WB |
| Immunogen | Recombinant protein of Human FIGF |
| Target | VEGFD |
| Preservatives | Store at -20°C or -80°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3. |
| Form/Appearance | liquid |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | 043915 |
| Dilution Range | WB: 1: 500-1: 2000, IHC-P: 1: 50-1: 200 |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of 293T cell lysat...