

## Product Datasheet

### Rad17 (phospho-Ser645) antibody (orb6829)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | Rad17 (phospho-Ser645) antibody   |
| <b>Species/Host</b>        | Rabbit  |
| <b>Conjugation</b>         | Unconjugated  |
| <b>Tested Applications</b> | ELISA, IF, IHC-Fr, IHC-P  |
| <b>Immunogen</b>           | KLH conjugated Synthesised phosphopeptide derived from human Rad17 around the phosphorylation site of Ser645 (PL(P-S)Q) |
| <b>Target</b>              | Phospho-Rad17 (Ser645)  |
| <b>Preservatives</b>       | 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.   |
| <b>Form/Appearance</b>     | Liquid  |
| <b>Concentration</b>       | 1mg/ml  |
| <b>Storage</b>             | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.  |
| <b>Note</b>                | For research use only   |
| <b>Isotype</b>             | IgG   |
| <b>Clonality</b>           | Polyclonal  |
| <b>MW</b>                  | 77 kDa  |
| <b>Uniprot ID</b>          | <a href="#">O75943</a>  |
| <b>Expiration Date</b>     | 12 months from date of receipt.   |