


Product Datasheet

Hepcidin-25 antibody (orb682923)

Description	Hepcidin-25 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant protein of Human HAMP
Target	HAMP
Preservatives	Store at -20°C or -80°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P81172
Dilution Range	WB: 1: 500-1: 2000, IHC-P: 1: 50-1: 200
Expiration Date	12 months from date of receipt.



Western blot analysis of MB435 cell line...