


Product Datasheet

MMP1 antibody (orb683093)

Description	MMP1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	A synthesized peptide derived from Human MMP1
Target	MMP1
Preservatives	0.01M PBS, pH 7.4 with 10 mg/ml rAlbumin and 0.1% Sodium azide
Form/Appearance	liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P03956
Dilution Range	WB: 1: 500-1: 2000
Expiration Date	12 months from date of receipt.



Western blot analysis of using MMP1 anti...