

Product Datasheet

ORF4 antibody (orb684616)

Description	ORF4 antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Peptide sequence around aa.9~13 (D-A-R-E-Q) derived from Cucumber necrosis virus RNA silencing suppressor.
Target	ORF4
Preservatives	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P15184
Dilution Range	WB: 1: 500-1: 3000
Expiration Date	12 months from date of receipt.