



Product Datasheet

FAM129A antibody (orb684749)

Description	FAM129A antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	peptide sequence around aa. 912~916 (K-E-G-E-G) derived from Human Niban.
Target	FAM129A
Preservatives	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9BZQ8
Dilution Range	WB: 1: 500-1: 1000, IF: 1: 100-1: 200
Expiration Date	12 months from date of receipt.



Western blot analysis of A431 and DU145 ...



Immunofluorescence analysis of HeLa cell...