



Product Datasheet

ADORA3 antibody (orb685275)

Description	ADORA3 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Fusion protein corresponding to a region derived from internal residues of Human Adenosine receptor A3
Target	ADORA3
Preservatives	Rabbit IgG in pH7.3 PBS, 0.05% NaN ₃ , 50% Glycerol.
Form/Appearance	liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P33765
Dilution Range	WB: 1: 500-1: 3000, IHC-P: 1: 50-1: 100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Human es...



Western blot analysis of Mouse testis ti...