


Product Datasheet

Histone H3K14me2 antibody (orb685610)

Description	Histone H3K14me2 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	A synthetic peptide of Human H3K14me2
Target	HIST3H3
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q16695
Dilution Range	WB: 1: 500-1: 2000, IHC-P: 1: 50-1: 100, IF: 1: 20-1: 50
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa and H3 pro...