

Product Datasheet

Defensin D2 antibody (Biotin) (orb685650)

Description	Defensin D2 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Plant
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Spinacia oleracea Defensin D2 protein (1-52AA)
Target	Defensin D2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P81571
Expiration Date	12 months from date of receipt.