

Product Datasheet

Allergen Ara h 1 antibody (HRP) (orb685652)

Description	Allergen Ara h 1 antibody (HRP)
Species/Host	Rabbit
Reactivity	Plant
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Peanut Allergen Ara h 1, clone P17 protein (26-216AA)
Target	Allergen Ara h 1, clone P17
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P43237
Expiration Date	12 months from date of receipt.