

Product Datasheet

GLU1 antibody (HRP) (orb685774)

Description	GLU1 antibody (HRP)
Species/Host	Rabbit
Reactivity	A. thaliana
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Arabidopsis thaliana Ferredoxin-dependent glutamate synthase 1, chloroplastic/mitochondrial protein (106-505AA)
Target	GLU1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9ZNZ7
Expiration Date	12 months from date of receipt.