

Product Datasheet

RB38 antibody (HRP) (orb685786)

Description	RB38 antibody (HRP)
Species/Host	Rabbit
Reactivity	Plant
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Chlamydomonas reinhardtii 38 kDa RNA-binding protein (1-382AA)
Target	RB38
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q6EMK7
Expiration Date	12 months from date of receipt.