

Product Datasheet

AOX1A antibody (Biotin) (orb685820)

Description	AOX1A antibody (Biotin)
Species/Host	Rabbit
Reactivity	A. thaliana
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Arabidopsis thaliana Ubiquinol oxidase 1a, mitochondrial protein (63-178AA)
Target	AOX1A
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q39219
Expiration Date	12 months from date of receipt.