


Product Datasheet

CESA4 antibody (orb685825)

Description	CESA4 antibody
Species/Host	Rabbit
Reactivity	A. thaliana
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Peptide sequence from Arabidopsis thaliana Cellulose synthase A catalytic subunit 4 [UDP-forming] protein (924-937AA)
Target	CESA4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q84JA6
Dilution Range	WB: 1: 100-1: 1000
Expiration Date	12 months from date of receipt.



Western blot analysis of Glycine max lea...