

Product Datasheet

CESA4 antibody (Biotin) (orb685828)

Description	CESA4 antibody (Biotin)
Species/Host	Rabbit
Reactivity	A. thaliana
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Peptide sequence from Arabidopsis thaliana Cellulose synthase A catalytic subunit 4 [UDP-forming] protein (924-937AA)
Target	CESA4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q84JA6
Expiration Date	12 months from date of receipt.