

Product Datasheet

B2m antibody (FITC) (orb686015)

| | |
|------------------------|--|
| Description | B2m antibody (FITC) |
| Species/Host | Rabbit |
| Reactivity | Mouse |
| Conjugation | FITC |
| Immunogen | Recombinant Mouse Beta-2-microglobulin protein (21-119AA) |
| Target | B2m |
| Preservatives | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol |
| Form/Appearance | Liquid |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | P01887 |
| Expiration Date | 12 months from date of receipt. |