



## Product Datasheet

### PPP5C antibody (orb686108)


<b>Description</b>	PPP5C antibody
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, IF, IHC, WB
<b>Immunogen</b>	Recombinant Human Serine/threonine-protein phosphatase 5 protein (370-488AA)
<b>Target</b>	PPP5C
<b>Preservatives</b>	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Uniprot ID</b>	<a href="#">P53041</a>
<b>Dilution Range</b>	WB: 1: 500-1: 2000, IHC-P: 1: 200-1: 500, IF: 1: 50-1: 200
<b>Expiration Date</b>	12 months from date of receipt.



Immunofluorescence analysis of A549 cell...



Western blot analysis of U87 whole cell ...



Immunohistochemical staining of human li...