

Product Datasheet

AOX1 antibody (Biotin) (orb686187)

Description	AOX1 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Plant
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Peptide sequence from Chlamydomonas reinhardtii Ubiquinol oxidase protein (202-217AA)
Target	AOX1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O65000
Expiration Date	12 months from date of receipt.