

Product Datasheet

KRT20 antibody (Biotin) (orb686695)

Description	KRT20 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Keratin, type I cytoskeletal 20 protein (207-424AA)
Target	KRT20
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P35900
Expiration Date	12 months from date of receipt.