

Product Datasheet

nef antibody (HRP) (orb686801)

Description	nef antibody (HRP)
Species/Host	Rabbit
Reactivity	Virus
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human immunodeficiency virus 1 Protein Nef protein (1-209AA)
Target	nef
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	A0A1L4CSJ1
Expiration Date	12 months from date of receipt.