



Product Datasheet

CACNA1B antibody (orb687036)

Description	CACNA1B antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC
Immunogen	Recombinant Human Voltage-dependent N-type calcium channel subunit alpha-1B protein (2132-2339AA)
Target	CACNA1B
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q00975
Dilution Range	IHC-P: 1: 20-1: 200, IF: 1: 50-1: 200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of SH-SY5Y C...



Immunohistochemical staining of human ki...