



Product Datasheet

GPR37 antibody (orb687084)


Description	GPR37 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Peptide sequence from Human Prosaposin receptor GPR37 protein (114-133AA)
Target	GPR37
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O15354
Dilution Range	WB: 1: 1000-1: 5000, IHC-P: 1: 500-1: 1000, IF: 1: 200-1: 500
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of PC-3 cell...



Western blot analysis of 293 whole cell ...



Immunohistochemical staining of human li...