



Product Datasheet

IL13RA1 antibody (orb687096)

Description	IL13RA1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC
Immunogen	Recombinant Human Interleukin-13 receptor subunit alpha-1 protein (22-234AA)
Target	IL13RA1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P78552
Dilution Range	IHC-P: 1: 20-1: 200, IF: 1: 50-1: 200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of U251 cell...



Immunohistochemical staining of human ki...