



Product Datasheet

Substance P antibody (orb687156)

Description	Substance P antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Protachykinin-1 protein (20-129AA)
Target	TAC1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P20366
Dilution Range	WB: 1: 1000-1: 5000, IHC-P: 1: 1000-1: 2000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human en...



Western blot analysis of HepG2 whole cel...