

Product Datasheet

hns antibody (FITC) (orb687390)

Description	hns antibody (FITC)
Species/Host	Rabbit
Reactivity	E. coli
Conjugation	FITC
Immunogen	Recombinant Escherichia coli DNA-binding protein H-NS protein (39-137AA)
Target	hns
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P0ACF8
Expiration Date	12 months from date of receipt.