

Product Datasheet

CLCN5 antibody (FITC) (orb687410)

Description	CLCN5 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human H(+)/Cl(-) exchange transporter 5 protein (553-746AA)
Target	CLCN5
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P51795
Expiration Date	12 months from date of receipt.