

Product Datasheet

VP2 antibody (HRP) (orb687521)

Description	VP2 antibody (HRP)
Species/Host	Rabbit
Reactivity	Virus
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Chicken anemia virus Dual specificity protein phosphatase VP2 protein (1-216AA)
Target	VP2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P69484
Expiration Date	12 months from date of receipt.