

## Product Datasheet

### Nucleoprotein antibody (Biotin) (orb687527)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | Nucleoprotein antibody (Biotin)  |
| <b>Species/Host</b>        | Rabbit   |
| <b>Reactivity</b>          | Virus  |
| <b>Conjugation</b>         | Biotin   |
| <b>Tested Applications</b> | ELISA  |
| <b>Immunogen</b>           | Recombinant Bovine respiratory syncytial virus Nucleoprotein protein (1-391AA) |
| <b>Target</b>              | N  |
| <b>Preservatives</b>       | Preservative: 0.03% Proclin 300.<br>Constituents: 50% Glycerol                 |
| <b>Form/Appearance</b>     | Liquid   |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C.<br>Avoid repeated freeze.               |
| <b>Note</b>                | For research use only  |
| <b>Isotype</b>             | IgG  |
| <b>Clonality</b>           | Polyclonal   |
| <b>Uniprot ID</b>          | <a href="#">P35943</a>   |
| <b>Expiration Date</b>     | 12 months from date of receipt.  |