

Product Datasheet

Hmgb1 antibody (FITC) (orb688042)

Description	Hmgb1 antibody (FITC)
Species/Host	Rabbit
Reactivity	Rat
Conjugation	FITC
Immunogen	Peptide sequence from Rat Hmgb1 protein (2-18AA)
Target	Hmgb1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P63159
Expiration Date	12 months from date of receipt.