


Product Datasheet

PCDHGB7 antibody (orb688220)

Description	PCDHGB7 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	Recombinant Human Protocadherin gamma-B7 protein (719-929AA)
Target	PCDHGB7
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9Y5F8
Dilution Range	IF: 1: 100-1: 500
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of SH-SY5Y C...