

Product Datasheet

RGS9BP antibody (FITC) (orb688266)

Description	RGS9BP antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Regulator of G-protein signaling 9-binding protein (1-210AA)
Target	RGS9BP
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q6ZS82
Expiration Date	12 months from date of receipt.