



Product Datasheet

Di-methyl-Histone H3(K27) antibody (orb688416)

Description	Di-methyl-Histone H3(K27) antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IP, WB
Immunogen	Synthetic Peptide
Target	Di-methyl-Histone H3(K27)
Preservatives	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Form/Appearance	liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P68431 , P84243 , Q71DI3
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa using Di-m...



Immunoprecipitation analysis of HeLa Cel...