


## Product Datasheet

### Di-methyl-Histone H3(K27) antibody (orb688417)

<b>Description</b>	Di-methyl-Histone H3(K27) antibody
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, WB
<b>Immunogen</b>	Synthetic Peptide
<b>Target</b>	Di-methyl-Histone H3(K27)
<b>Preservatives</b>	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
<b>Form/Appearance</b>	liquid
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal
<b>Uniprot ID</b>	<a href="#">P68431</a> , <a href="#">P84243</a> , <a href="#">Q71DI3</a>
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of 1) HeLa, 2) Rat...