

Product Datasheet

Tri-methyl-Histone H3(K4) antibody (orb688418)

Description	Tri-methyl-Histone H3(K4) antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthetic Peptide
Target	Tri-methyl-Histone H3(K4)
Preservatives	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Form/Appearance	liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P68431, P84243, Q71DI3
Expiration Date	12 months from date of receipt.



Western blot analysis of 1) HeLa, 2) Rat...