



Product Datasheet

Phosphoserine antibody (orb688486)


Description	Phosphoserine antibody
Species/Host	Mouse
Reactivity	All
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Purified Protein
Target	Phosphoserine
Preservatives	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Form/Appearance	liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Dilution Range	WB:1:500-1:2000, IHC-P:1:100-1:300
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Human Br...



Western blot analysis of 1) Hela, 2)Rat ...



Immunohistochemical staining of human lu...