

## Product Datasheet

### TUBB antibody (HRP) (orb688489)

<b>Description</b>	TUBB antibody (HRP)
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Canine, Gallus, Hamster, Human, Insect, Monkey, Mouse, Rabbit, Rat, Sheep, Yeast
<b>Conjugation</b>	HRP
<b>Tested Applications</b>	ELISA, IHC, WB
<b>Immunogen</b>	Recombinant Protein
<b>Target</b>	TUBB
<b>Preservatives</b>	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
<b>Form/Appearance</b>	liquid
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal
<b>Uniprot ID</b>	<a href="#">P07437</a>
<b>Dilution Range</b>	WB:1:500-1:5000, IHC-P:1:50-1:500
<b>Expiration Date</b>	12 months from date of receipt.