

## Product Datasheet

### CD4/CD8 antibody (FITC/PE) (orb688640)

<b>Description</b>	CD4/CD8 antibody (FITC/PE)
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, FC, IF
<b>Immunogen</b>	Recombinant Protein
<b>Target</b>	CD4/CD8
<b>Preservatives</b>	Phosphate-buffered solution, pH 7.4, containing 0.09% sodium azide and 0.2% (w/v) rAlbumin
<b>Form/Appearance</b>	liquid
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG2b/IgG1
<b>Clonality</b>	Monoclonal
<b>Uniprot ID</b>	<b>P01730</b>
<b>Expiration Date</b>	12 months from date of receipt.