

Product Datasheet

CD4/CD8/CD3 antibody (FITC/PE/Cy5) (orb688641)

Description	CD4/CD8/CD3 antibody (FITC/PE/Cy5)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF
Immunogen	Recombinant Protein
Target	CD4/CD8/CD3
Preservatives	Phosphate-buffered solution, pH 7.4, containing 0.09% sodium azide and 0.2% (w/v) rAlbumin
Form/Appearance	liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG2b/IgG1/IgG2a
Clonality	Monoclonal
Uniprot ID	P01730
Expiration Date	12 months from date of receipt.