

## Product Datasheet

### PDHA1 antibody (orb688688)

<b>Description</b>	PDHA1 antibody
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Canine, Human, Mouse, Porcine, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, ICC, WB
<b>Immunogen</b>	Recombinant fragment , corresponding to amino acids 131aa-389aa of human PDHA1 (C-terminus).
<b>Target</b>	PDHA1
<b>Preservatives</b>	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol
<b>Form/Appearance</b>	liquid
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Uniprot ID</b>	<a href="#">P08559</a>
<b>Expiration Date</b>	12 months from date of receipt.