


## Product Datasheet

### Rubisco (Large Chain) antibody (orb688722)

<b>Description</b>	Rubisco (Large Chain) antibody
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Plant
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, WB
<b>Immunogen</b>	Protein
<b>Target</b>	Rubisco (Large Chain)
<b>Preservatives</b>	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of Arabidopsis wit...