



Product Datasheet

TUBE1 antibody (orb688741)

Description	TUBE1 antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Synthetic Peptide of Epsilon Tubulin
Target	TUBE1
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	Q9UJT0
Dilution Range	IHC-P:1:100-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Human Co...



Immunohistochemical staining of human CO...