



Product Datasheet

CD146 antibody (orb688859)


Description	CD146 antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, WB
Immunogen	Recombinant Human Cell surface glycoprotein MUC18 protein (50-646AA)
Target	CD146
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	P43121
Dilution Range	WB:1:500-1:2000, IF:1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HeLa cell...



Immunofluorescence analysis of MCF-7 cel...



Immunofluorescence analysis of THP-1 cel...