



Product Datasheet

CD19 antibody (orb688865)


Description	CD19 antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IHC, WB
Immunogen	Recombinant Human CD19 protein (20-291AA)
Target	CD19
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P15391
Dilution Range	WB:1:2000-1:80000, IHC-P:1:50-1:200, IF:1:50-1:200
Expiration Date	12 months from date of receipt.



Flow cytometric analysis of Raji cells u...



Western blot analysis of Rat Heart Tissue...



Western blot analysis of Raji whole cell...