



Product Datasheet

Beta galactosidase antibody (orb688869)


Description	Beta galactosidase antibody
Species/Host	Mouse
Reactivity	E. coli
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IP, WB
Immunogen	Recombinant Escherichia coli Beta-galactosidase protein
Target	lacZ
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P00722
Dilution Range	WB:1:1000-1:64000, IF:1:50-1:200, IP:1ul-4ul, FC:1:100-1:300
Expiration Date	12 months from date of receipt.



Flow cytometric analysis of 293T transfe...



Western blot analysis of Human Colon Car...



Immunofluorescence analysis of 293T cell...