


## Product Datasheet

### CD9 antibody (orb688874)

<b>Description</b>	CD9 antibody
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, IHC, WB
<b>Immunogen</b>	Recombinant Human CD9 antigen protein (112-195AA)
<b>Target</b>	CD9
<b>Preservatives</b>	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG2a
<b>Clonality</b>	Monoclonal
<b>Uniprot ID</b>	<a href="#">P21926</a>
<b>Dilution Range</b>	WB: 1:1000-1:5000, IHC-P:1:500-1:1000, IF:1:150-1:300, FC:1:50-1:100, IP:1 µl-4 µl
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of Rat Heart Tissue...