



Product Datasheet

YWHAZ antibody (orb688879)


Description	YWHAZ antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IHC, WB
Immunogen	Recombinant Human 14-3-3 protein zeta/delta protein (133-212AA)
Target	YWHAZ
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Uniprot ID	P63104
Dilution Range	WB:1:5000-160000, IHC-P:1:50-1:200, IF:1:50-1:200, FC:1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of A549 cell...



Immunofluorescence analysis of HeLa cell...



Immunofluorescence analysis of U251 cell...