



Product Datasheet

HA-Tag antibody (orb688917)


Description	HA-Tag antibody
Species/Host	Mouse
Reactivity	All
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IP, WB
Immunogen	YPYDVPDYA synthetic peptide conjugate to KLH
Target	HA-Tag
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Dilution Range	WB:1:5000-1:160000, IF:1:50-1:200, IP: 1:5000-1:100000, FC:1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of fusion protein ...



Western blot analysis of Human Liver Car...



Western blot analysis of using HA-Tag an...