

## Product Datasheet

### Goat Anti-Mouse IgG(H+L) antibody (Biotin) (orb688926)

<b>Description</b>	Goat Anti-Mouse IgG(H+L) antibody (Biotin)
<b>Species/Host</b>	Goat
<b>Reactivity</b>	Mouse
<b>Conjugation</b>	Biotin
<b>Tested Applications</b>	ELISA, FC, IF, IHC, IP, WB
<b>Immunogen</b>	Mouse IgG(H+L)
<b>Preservatives</b>	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Dilution Range</b>	WB:1:5000-20000, IHC-P:1:200-2000, IF:1:200-2000, IP:1:200-1:2000, FCM:1:200-2000, ELISA: 1:5000-20000
<b>Expiration Date</b>	12 months from date of receipt.